User Manual

Ultra-Stable Temperature
Chamber (Air Bath)

TK-105 US

SAMPLE

Kambič

WARRANTY

Kambic d.o.o. guaranties that the product delivered has been thoroughly tested to ensure that it meets its published specifications. The warranty is valid only if the product has been installed and used in accordance to the instructions supplied by Kambic d.o.o.

Kambic d.o.o. shall in no event be liable for incidental or consequential damages, including without limitations, lost profits, loss of income, loss of business opportunities, loss of use, and other related exposures, caused by incorrect use of product.

Warranty is valid for **13 months** from date of shipping and includes spare parts and labor.

Type:	TK-105 US
Serial number:	
Date of shipping:	
Signature/stamp:	

Manufacturer:

Kambic d.o.o.
Metliska cesta 16
8333 Semic
Slovenia

T +386 7 35 65 220

F +386 7 35 65 232

E info@kambicmetrology.com

W www.kambicmetrology.com

Thank you for your trust and for buying this device. We hope it will successfully serve the purpose for many years.

**CAREFULLY READ THESE INSTRUCTIONS
BEFORE USING THE DEVICE**

SAMPLE

CONTENT

1. WARNINGS AND SAFETY PRECAUTIONS	4
2. INTRODUCTIONS	14
3. TECHNICAL DATA	15
4. INSTRUCTIONS FOR USE	16
4. 1 Basic Setup	16
4.2 Control Panel With Touch Display	16
4.2.1 Main Screen	17
4.2.2 Info Screen	17
4.2.3 Unlock Screen	18
4.2.4 Graph Screen	18
4.2.5 Main Settings	19
4.2.6 Chamber-Pid Settings	19
4.2.7 Thermoelectric Cooler-Pid Settings	20
4.2.8 Temperature Corrections.....	20
4.2.9 Stop Screen	21
4.2.10 Door Alarm Screen.....	21
5. SERVICE AND MAINTENANCE	22
6. REQUIRES SPACE AND CONECTIONS	24
7. DISMANTLING AND DISPOSAL	25
8. ABBRIVATIONS	26
9. TROUBLESHOOTING	27
10. RS 232 communication protocol	28
11. WIRING DIAGRAM	
12. DECLARATIONS	

EN

1. WARNINGS AND SAFETY PRECAUTIONS

Always keep this manual in close reach for reference. This product should only be used and operated by skilled and trained personnel familiar with technical capabilities and potential safety hazards.

Always follow and obey the national regulations related to technology and labor laws.

To avoid injuries and equipment damage study and follow the safety instructions and warnings.

Legal:

The content of the manual is subject to change without notice.

Understanding and following this operating manual is first step to hazard-free and safe use of the equipment.

This user manual can not reach out to all application that might even develop in the future. If you find information unclear or you are not sure about safe use of the product in your specific case, do not hesitate to contact your dealer or us via info@kambic.com for further information or opinions.

In no event shall Kambic d.o.o. be held liable for any damages, direct or incidental related to the use of this product and manual.

Product can be equipped with variety of accessories and options. This warning list might include elements that are applicable only with certain models and accessories.

READ ALL INSTRUCTIONS BEFORE USE!

Reduce the possibility of injury and failure.

WARNING

Before you perform any kind of maintenance, disconnect all electrical power to the unit! DO NOT operate this unite without a properly grounded, properly polarized power cord. DO NOT connect this unite to a non-grounded, non-polarized outlet.

DANGER

Appliance must not be installed in a room where the possibility of explosive atmosphere is present!

DANGER OF BREAKDOWN

Service only allowed to authorized personnel!

DANGER OF FIRE

Heating of flammable materials is not allowed!

WARNING

Close the doors tightly before start! Do not open doors during the cycle!

DANGER OF BURNS

Use safety gloves when operating with high/low temperatures!

EXPLOSION HAZARD

Do not process any materials that could generate explosive hazardous Atmosphere!

EXPLOSION HAZARD

Do not operate the equipment in explosive hazardous atmosphere!

ELECTRIC HAZARD

Do not operate the equipment in wet and condensing environment!

STABILITY HAZARD

Do not load or affect product stability with heavy loads or climbing the product!

TRANSPORT HAZARD

Transport the chamber only in original packing!
Always secure the chamber stability during reposition!
Never use doors or handles or accessories for holding point during reposition!

OVERHEATING HAZARD

Do not install the equipment in non ventilated environments!
Ensure appropriate space and volume around the chamber!
Ensure sufficient heat dissipation for product cooling!

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2. INTRODUCTIONS

Ultra-stable temperature chamber (air bath) is design to maintain the uniform constant temperature over a wide operating range.

Air Baths are a superb choice for precision controlled operating environments without the mess and expense of fluid bath operation.



The device housing is made of powder coated aluminized sheet metal. Internal chamber is made of stainless steel AISI 304. Inside are two shelves witch allowing easy access to the inserted objects. Easy access to the inside of the chamber is achieved through the conveniently mounted front door. Wires or cables are brought into the chamber through a small, insulated opening on the side of the unit.

The integrated control unit adds convenience as it eliminates a separate controller. The internal temperature is shown on an external display located on the front of the chamber. A touch display can be used to set the individual temperatures between 15 °C and 50 °C.

Cooling and heating for the chamber is provided by a Peltier Cooling System. Internal fan is used to circulate the air inside the chamber to provide temperature stability and uniformity.

3. TECHNICAL DATA

TK-105 US	
External dimensions	705 x 1130 x 700mm (W x H x D)
Internal dimensions	400 x 640 x 410mm (W x H x D)
Volume	105 L
Temperature range	+15°C ...+50°C
Temperature display resolution	0.001°C
Temperature set resolution	0.001°C
Temperature stability	+/- 0.005°C @ 23°C
Control	MPC-PID
Power supply	230V±10%; 50Hz, 60Hz
Wattage	200W
Interface	RS-232
Shelve	2 included
Access port	2 x 55x20mm included
Weight	112kg

4. INSTRUCTIONS FOR USE

Step 1: Connect the power cord into the power supply 230V±10%; 50Hz, 60Hz

Step 2: Close the chamber door.

Step 3: Put the main switch in position 1 (on the back of the device)


Step 3: Press On/off switch

Step 4: Select the desired temperature and press "**START button**"


Step 5: From this moment on the device runs completely automatically

Step 6: Press "**STOP**" to stop the device

4.2 Control panel with touch display


4.2.1 Main screen


4.2.2 INFO screen

Press INFO button to see this screen:

4.2.3 Unlock screen

Press Unlock button to see this screen:

Enter the password 0207 and confirm with OK button. That will enable editing PID settings and temperature corrections.

To disable the possibility of editing PID setting and temperature corrections press the "Lock" button. The controller will lock itself automatically after 2 minutes of no action.

4.2.4 Graph screen

Press Graph button to see this screen:

4.2.5 Main settings

The screenshot shows the 'MAIN SETTINGS' screen with the following items:

- Temperature:** 35.00°C (with a red hand icon pointing to the value)
- Vent.SpeedCtrl:** ON
- LCD Backlight:** 5
- Four directional arrow buttons (Left, Down, Up, Right) at the bottom.

Callouts from the right side of the image point to the following elements:

- Touch value to select:** Points to the temperature value.
- Set temperature:** Points to the temperature value.
- Fan speed control:** Points to the 'ON' status. Below it, text reads: 'OFF: fan operates on full RPM', 'ON: fan RPMs are controlled'.
- Set LCD brightness:** Points to the '5' value.
- Next screen:** Points to the Right arrow button.
- Increase value:** Points to the Up arrow button.
- Decrease value:** Points to the Down arrow button.
- Previous (main) screen:** Points to the Left arrow button.

4.2.6 Chamber-PID settings

The screenshot shows the 'Chamber-PID SETTINGS' screen with the following items:

- Kp-Heating Gain:** 2.0 (with a red hand icon pointing to the value)
- Kp1-Cooling Gain:** 2.0
- Ti-Integral time:** 150.0
- Td-Derivat.const.:** 0.0
- (Fact.Defaults: Kp= 2.0,Kp1= 2.0,Ti=150,Td= 0)
- Four directional arrow buttons (Left, Down, Up, Right) at the bottom.

Callouts from the right side of the image point to the following elements:

- Touch value to select:** Points to the Kp-Heating Gain value.
- Next screen:** Points to the Right arrow button.
- Increase value:** Points to the Up arrow button.
- Decrease value:** Points to the Down arrow button.
- Previous screen:** Points to the Left arrow button.

Default values	
Kp-Heating Gain	2.0
Kp1-Cooling Gain	2.0
Ti-Integral time	150.0
Td-Derivat.const.	0.0

4.2.7 Thermoelectric cooler - PID settings

TEC - PID SETTINGS

Kp-Heating Gain: 30.0
 Kp1-Cooling Gain: 60.0
 Ti-Integral time: 50.0
 Td-Derivat.const.: 0.0
 (Fact.Defaults: Kp= 30.0,Kp1= 60.0,Ti= 50,Td= 0)

Navigation buttons: Left, Down, Up, Right

Labels: Touch value to select, Next screen, Increase value, Decrease value, Previous screen

Default values	
Kp-Heating Gain	30.0
Kp1-Cooling Gain	60.0
Ti-Integral time	50.0
Td-Derivat.const.	0.0

4.2.8 Temperature corrections

Temp.Corr. T= 34.973°C

Corr. at 0.0°C = 0.000
 Corr. at 15.0°C = 15.000
 Corr. at 25.0°C = 25.000
 Corr. at 40.0°C = 40.000


Navigation buttons: Left, Down, Up, Right

Labels: Touch value to select, Next screen, Increase value, Decrease value, Previous (main) screen

Default values	
Corr. At 0.0°C	0.000
Corr. At 15.0°C	14.980
Corr. At 25.0°C	24.870
Corr. At 40.0°C	39.980

4.2.9 STOP screen

If you press stop button (1 second or more) to stop the cycle you will see this screen:

4.2.10 Door alarm screen

If you open the door while cycle is in progress there will be an audio alarm and you will see this screen:

5. SERVICE AND MAINTENANCE

WARNING

In case of maintenance interventions, disconnect power supply!

Maintenance and service within the warranty period may only be performed by the authorized service person.


Please advise Kambic in before alternation & modification.

All spare parts are in stock in the producer's magazine. It is only allowed to use the original spare parts.


Keep the inside and outside surfaces of the chamber clean. For cleaning use some soft cleaning detergents and soft cloth.

Clean the rear cooling fan and cooling ribs.

Check and remove condensate under rear cooling ribs.

Check and remove condensate from internal chamber.

internal chamber opening

In the event of failure, call the service!

Maintenance plan

OPERATION	FREQUENCY
Clean the rear cooling fan and cooling ribs	monthly
Check and remove condensate from internal chamber	weekly
Check and remove condensate under rear cooling ribs	weekly
General cleaning of the apparatus housing	weekly

6. REQUIRES SPACE AND CONNECTIONS


The appliance must be positioned at least 300 mm away from a wall to allow the cooling of Peltier cooling system via its rear fan.

Additional space or relocation of the apparatus may be necessary when carrying out any maintenance or servicing.

The user must provide electrical connection: 200W / 230V \pm 10% / 50Hz, 60Hz.

Environmental and operational conditions

Temperature	+23°C \pm 5°C
Relative humidity	from 15 % to 75 %

7. DISMANTLING AND DISPOSAL

This symbol means do not dispose of your product with your other household waste. Instead, you should protect human health and the environment by handing over your waste equipment to a designated collection point for the recycling of waste electrical and electronic equipment. For more information please contact your household waste disposal service.

8. ABBRIVATIONS

T	chamber temperature
Tset	set temperature
Ttec	thermoelectric cooler temperature
TecPwr	thermoelectric cooler power
PID	proportional–integral–derivative (control loop)
LCD	Liquid Cristal Display
LED	light-emitting diode
RS-232	standard communication interface

9. TROUBLESHOOTING

Alarm door opened

- ➔ check the end switch
- ➔ check the locking mechanism

Water drip from chamber

- ➔ remove condensate
- ➔ check the chamber seal
- ➔ check the access port seal

Low temperature stability & uniformity

- ➔ check the door seal
- ➔ check the operation of the rear fan
- ➔ check the operation of the inner fan
- ➔ check the cooling ribs inside & outside
- ➔ check the temperature sensor

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- ➔ replace it with the original

Damaged foam on access port

- ➔ replace with a closed cell type foam

Controller malfunction

- ➔ call service

10. RS-232 COMMUNICATION PROTOCOL

COM port settings:

- 9600 Baud
- No parity
- 8 Data bit
- 1 Stop bit

Connector pins info:

- 3 = TXD RS232 output from controller
- 7 = GND
- 11 = RXD RS232 input to controller

All commands must be written in capital letters and identical as written in this description!

After the receipt of the command (START, STOP, SETTEMP) the data is immediately written in the eeprom. This way the controller will always return to the last settings in case of a power failure.

No.	Command:	Action:
1.	STOP <LF>	Operation of the bath stopped. Equal to STOP on the controller.
2.	START<LF>	Operation of the bath activated. Equal to the START on the controller.
3.	SETTEMP= 37.0<LF>	Set temperature with single decimal
4.	SETTEMP?<LF>	Request of set temperature. Will return current set temperature. SETTEMP=37.000<CF><LF>
5.	TXD<LF>	Request of set and current temperature. Will return current set temperature and current measured temperature. Tm=37.001 Ts=37.0<CR><LF>
6.	KP=120<LF>	Set P component in PID controller
7.	KP?<LF>	Request of current P component. Will return KP=120<CR><LF>
8.	TI=600<LF>	Set I component in PID controller
9.	TI?<LF>	Request of current I component. Will return TI=600<CR><LF>
10.	TD=0<LF>	Set D component in PID controller
11.	TD?<LF>	Request of current D component. Will return TD=0<CF><LF>
12.	TRIG=10<LF>	Automatic RS232 output period 10 seconds. 0 seconds automatic output disabled. Will output: 10.06.3 Tm=37.001 Ts=37.0 P=75WD=0<CR><LF> Usable with combination to dot printer EPSON LX-300
13.	TRIG?<LF>	Request of current output period. Will return TRIG=10<CR><LF>